

# BUEHLER® SurfMet® 1 METALLOGRAPHIC AND SPECTROGRAPHIC BELT SURFACERS

- Recessed 4" x 9" (10 x 23cm) Grinding Tables with Belt Guards
- Built-in Flushing Systems and Stainless Steel Sinks
- Fast, Easy Belt Installation with Non-slip Belt Mechanism and Fingertip Centering
- Corrosion-resistant Stainless Steel Belt Compartments and Piping
- Convenient On/Off Switch with Indicator Light
- Accessory Recirculating System for Wet Grinding
- Accessory Switch Activates Outlet on Back
- Water Flow Control Knob



**Buehler SurfMet® I Belt Surfacers** make grinding metallographic or spectrographic samples easy and efficient. Features 1/3 horsepower (249W), heavy-duty, fused motors for superior durability and performance. Replaceable stainless steel wear plates are designed to maximize machine life. Separate belt flushing and coolant controls guarantee adequate removal of particulate matter from the grinding belt, for efficient cutting action, while providing critical specimen cooling. The stainless steel sink with faucet allows the convenience of quickly grinding and washing the specimen in one location; this feature is particularly useful between coarse and fine grinding steps and for production applications.

The SurfMet I Belt Surfacers are offered as tabletop models mounted on sturdy Tech-met® tables with chemical-resistant laminate tops. The optional recirculating system utilizes a corrosion inhibiting lubricant that extends machine life, provides faster, cleaner grinding, and collects swarf for environmentally correct disposal.

## Specifications

**16-1270-160 SurfMet I Belt Surfacers**, corrosion-resistant modified epoxy coated cabinet, 4" x 9" (10.2cm x 22.9cm) grinding area, cam operated belt centering device, belt coolant, and compartment flushing system, stainless steel sink with faucet, illuminated on-off switch, 115 V/60Hz/1 phase, 1/3 Hp (249W) fused motor. Includes three 16-5200 CarbiMet MetSplice Abrasive Belts and operating instructions. Dimensions: 18" W x 21" D x 11" H (46cm x 53cm x 28cm). Shipping weight: 95 lbs. (43kg).

**16-1270-250 SurfMet I Belt Surfacers**, same as 16-1270-160, but for operation on 220-240V/50Hz/1 phase.

## Accessories

**10-1250-160 Recirculating System**, including baffled and compartmented tank, 115V/50-60Hz/1 phase, 1/16 Hp (50W) motor connected directly to pump. Dimensions: 13" W x 21" D x 16" H (3cm x 53cm x 41cm). Shipping weight: 75 lbs. (34kg).

**10-1250-250 Recirculating System**, same as 10-1250-160, but for operation on 220-240V/50Hz/1 phase.

**16-1279 Copper Wear Plate Retrofit Kit**, for Dry Grinding Applications.

**16-1278 Wear Plate Replacement Kit**, for Belt Surfacers.

**16-3170 Hose and Connector Assembly**, for connection of SurfMet I tabletop model to Recirculating System.

**16-1499 TechMet Cabinet**, for use with all Buehler tabletop belt surfacers. Includes heavy duty TechMet metal base cabinet, chemical-resistant laminate top, with plumbing and electrical service access holes, leveling feet for ease of installation. Overall dimensions: 28 1/2" W x 27 1/4" D x 36" H (72 x 70 x 91cm). Approximate shipping weight: 132 lbs. (60kg).

\* Contact Buehler for dry grinding operation recommendations.

For a complete listing of consumables, please refer to our Consumables Buyer's Guide. Buehler continuously makes product improvements; therefore, technical specifications are subject to change without notice.

© 2010 BUEHLER®, a division of Illinois Tool Works, Inc. Printed in U.S.A. 00M0410 FN00535 Rev.3 \*For metallurgical consumables produced by BUEHLER®



## BUEHLER

**BUEHLER®**, a division of Illinois Tool Works, Inc.  
– **Worldwide Headquarters**  
41 Waukegan Road • P.O. Box 1  
Lake Bluff, Illinois 60044-1699 USA  
Tel: (847) 295-6500 • Fax: (847) 295-7979  
Sales: 1-800-BUEHLER (1-800-283-4537)  
www.buehler.com  
Email: info@buehler.com

### BUEHLER GMBH - European and MESA Headquarters

In der Steele 2 • 40599 Düsseldorf  
Postfach 16 03 55 • 40566 Düsseldorf  
Telefon: (49) 211 974100 • Telefax: (49) 211 9741079  
www.buehler-met.de  
Email: info@buehler-met.de

### BUEHLER FRANCE

Téléphone: 0800 89 73 71  
Télécopie: 0800 88 05 27  
www.buehler.fr  
Email: info@buehler.fr

### BUEHLER UNITED KINGDOM

Telephone: 0800 707 6273  
Fax: 0800 707 6274  
www.buehler.co.uk  
Email: sales@buehler.co.uk

### BUEHLER CANADA

10 Carlow Court, Unit #2  
Whitby, Ontario L1N 9T7  
Telephone: (905) 430-4684  
Fax: (905) 430-4647  
Sales Telephone: 1-800-268-3593  
Email: info@buehler.ca

### BUEHLER, ASIA-PACIFIC

5/F Vogue Centre  
696 Castle Peak Road  
Lai Chi Kok, Kowloon  
Hong Kong, SAR, China  
Telephone: (852) 2307 0909  
Fax: (852) 2307 0233